NIRT Weekly Barley SRS Worksheet

Location: Instrument Constants:			Serial# "O"			_	Date:				
Protein Constants:			Slope			<u>—</u>	Intercept				
Date:			1				1		1		
Tech:							_				
Temp:											
Humidity:											
Trummanty.					PROTE	INI					
SRS	Baseline	Result	Diff	Result	Diff	Result	Diff	Result	Diff	Result	Diff
1											
1											
2											
2											
3											
3											
4											
4											
5											
5											
Average											
Range of Differences (Max.		x. 0.60)									
			1								
			2nd Previous Run								
Bias Calculation:			Previous Run								
One Run (.12)											
Two Runs (.09)											
Three Runs (.06)											
Five Runs (all +/-)(.04)											