NIRT Barley SRS Bias Log

Tecator Model:	Serial # :	Location:
NOTE: Only results less th	an two weeks old and within a 5 °F to	emperature range may be averaged. Results
obtained after a bias adjus	tment, SRS replacement, or NIRT repo	air or replacement cannot be average with results
		be applied if five consecutive, valid runs are all
positive or all negative). E	Every effort should be made to maintai	in the bias within the 0.06% Protein tolerance.

		Protein Class:								
			Today's Multi-Run Averages Bias 2-Run 3-Run 5-Run							
			Bias	2-Run	3-Run	5-Run	B*			
Date mm/dd/yy	Temp	RH	±0.12	±0.09	±0.06	±0.04	R*	Oper Init	Comments	
*(B) bias adju	*(B) bias adjustment, or (R) SRS replaced. If the NIRT is repaired or replaced, start a new log.									