

*NOTE: Only results less than two weeks old and within a 5°F temperature range may be averaged. Results obtained after a slope or bias adjustment, SRS replacement, or NIRT repair or replacement cannot be averaged with results obtained prior to the change (also, the 5-Run correction can only be applied if five consecutive, valid runs are all positive or all negative). Every effort should be made to maintain the bias within the 0.05% Protein tolerance.*

Date mm/dd/yy	Temp							Oper Init	Comments
		RH	Today's Bias	Multi-Run Averages			S*		
				2-Run	3-Run	5-Run	B*		
			$\pm 0.10$	$\pm 0.07$	$\pm 0.05$	$\pm 0.03$	R*		

\* (S) slope adjustment, (B) bias adjustment, or (R) SRS replaced. If the NIRT is repaired or replaced, start a new log.