Weekly NIRT Corn SRS Worksheet

Constituent: Date: Temperature:			_ Lo	ocation:	Serial #:						
		_ Operator:									
Instrument Constants: Slope:			"O":		"P":						
Date:											
Temp:	:										
SRS	Baseline	Result	Diff	Result	Diff	Result	Diff	Result	Diff	Result	Diff
1											
1											
2											
2											
3											
3											
4											
4											
Avg											
Range											

Transfer bias results to the SRS Bias Log, calculate multi-run averages, and apply tolerances. Adjust intercept constant for a constituent if it is out of tolerance.