

*NOTE: Only results less than two weeks old and within a  $\pm 5^{\circ}\text{F}$  temperature range may be averaged. Results obtained after a slope or bias adjustment, SRS replacement, while outside of the RH limit, or NIRT repair or replacement can not be averaged with results obtained prior to the change (also, the 5-Run correction can only be applied if five consecutive, valid runs are all positive or all negative). Every effort should be made to maintain the bias within the  $\pm 0.07$  percent protein,  $\pm 0.07$  percent oil, and  $\pm 0.20$  percent starch tolerances.*

[illegible]

\* (S) slope adjustment, (B) bias adjustment, or (R) SRS replaced. If the NIRT is repaired or replaced, start a new log.